

Helps the patient to have a bath or shower without getting their lower leg or foot dressing wet.



Prescription Request Form

Patient _____ Hospital / Clinic _____

Clinician _____ Date _____

Half Leg LimbO - Seals just above the knee, ideal for lower leg bandages, casts or bulky foot dressings. For bath or shower. Sizing Instructions: Check your patient's height and weight to select the correct size half leg LimbO.



Patient Height	Patient Weight	Model	PIP Code	✓
5'5" (165 cm) and above	Up to 10 Stone (64kg)	MP76 Slim Leg	333-7581	
	10 to 16 Stone (64 to 102kg)	MP80 Average Build	320-0045	
	16 to 22 Stone (102 to 140kg)	MP180 Large Leg	320-0060	
Under 5'5" (165cm)	Up to 10 Stone (64kg)	MP76S Slim Leg Short	333-7599	
	10 to 16 Stone (64 to 102kg)	MP80S Short Leg	320-0052	
	16 to 22 Stone (102 to 140kg)	MP180S Large Leg Short	320-0078	

Foot LimbO - Seals just above the ankle, ideal for light dressings on the foot. Shower use only.



Leg Circumference (5cm above ankle bone)	Model	Size	PIP Code	✓
20 - 25cm (8-10")	MP20	Small	416-6476	
23 - 34cm (9-13")	MP25	Medium to Large	416-6484	

Please note: Maximum UK shoe size is 11. If larger, consider a half leg LimbO option.

Advice line: **01243 578755** Email: info@limboproducts.co.uk www.limboproducts.co.uk

Thesis Technology Products Ltd., Brooks Green Farm, Brooks Lane, Bosham, West Sussex, PO18 8JX

LFPP0321