

The original waterproof protector

Helping patients shower and bathe with confidence

Waterproof Protection for Casts, Air Walker Boots, Dressings and PICC Lines/Midlines



# Help restore patients' well-being with a LimbO waterproof protector

Thank you for your professional interest in the LimbO Waterproof Protector. This guide explains why the LimbO is recognised as a leading waterproof protection solution for patients and quickly assists you with product selection. The guide highlights how to obtain a LimbO including options on FP10 prescription or the NHS Shared Business Services framework (NHS SBS). We have also included a section about our outdoor cast protectors, cast accessories, mobility aids and a range of

For further support and resources please visit www.limboproducts.co.uk o

## **Contents**

Benefits of the LimbO	3
Adult LimbOs	4 - 7
Elbow LimbOs (adult & child)	5
Child LimbOs	8 - 9
Cast / Dressing / PICC Line Accesso	ries 10
Outdoor Cast / Dressing Protectors	11
Daily Living Aids and Diabetic Socks	12
Body Block	13
Size Guide	14 - 15
Contact Details	Back cover

"Patient use of the LimbO reduces the chance of wet dressings and casts and the associated problems of discomfort, clinical complications and delays to recovery".

# Why is the LimbO the most recommended protector?

## Helping the Patient

- More independence during recovery
- Easy to use
- Maintain personal hygiene
- Durable and submergible
- Reliable waterproof protection
- Adult and child sizes
- Easy to order via phone, web or Amazon
- Prices from £9.15
- Rapid delivery
- Latex free
- Intended for single patient use

## **Helping Your Team**

- Reduce healthcare risks associated with wet dressings and wet casts
- Reduce avoidable NHS costs resulting from re-casting and re-dressing
- Easy size selection
- Certain models available on FP10 prescription and NHS SBS framework.
- Bespoke sizes available for unique cases
- 27 years as a trusted solution
- Dedicated practitioner content on our website

## Featuring the unique LimbO neoprene seal

- Flexible and comfortable
- Superior design providing a complete waterproof barrier against the skin
- Gentle non-constrictive seal, made with premium neoprene
- Guaranteed for twelve weeks with regular use
- Proudly designed and manufactured in the UK



# Guaranteed waterproof protection for casts, dressings or PICC lines



## Adult Half Arm

For casts or dressings below the elbow.

- Seals just above the elbow
- Suitable for the bath or shower
- NHS SBS framework options



Height (cm)	Weight (kg)	Length (cm)	Circumfer Limb	rence (cm) Cast	Model	NHS SBS
152+	Under 51	49	22-25	33	M50A	
152+	51-70	53	25-29	40	M60	✓
152+	70+	53	30-39	44	M67	✓



## Adult Full Arm

For casts or dressings covering the whole arm.

- Seals at the top of the upper arm
- Suitable for the bath or shower
- NHS SBS framework options



	Height (cm)	Weight (kg)	Length (cm)	Circumference (cm) Limb Cast		Model	NHS SBS
Ī	152+	51-70	69	25-29	40	M60L	✓
	152+	70-108	69	30-39	44	M70	✓
	152+	108+	69	39-54	57	M87	✓




## Adult and Child Elbow

For PICC Lines, Midlines, or dressings on the elbow or upper arm.



- Ideal to support patients duringIV treatment
- Suitable for the bath or shower
- Prescription options
- Child option



39-54



152+

## **Adult Hand**

108+

For light dressings on the hand.

- Seals on the forearm
- Not suitable for plaster casts or following thumb surgery. Select Adult Half Arm model instead



	Height (cm)	Weight (kg)	Length (cm)	Circumf Limb	ference (cm) Cast	Model
Ī	Slim to Medium Build		30	22-25	Bandage Only	MITT

# The comfortable, trusted solution for staying dry in the shower or bath



## Adult Half Leg

For casts or dressings below the knee.

- Seals just above the knee
- Suitable for use with diabetic ulcers
- Suitable for the bath or shower
- Prescription options
- NHS SBS framework options



Height (cm)		Length (cm)					NHS SBS
165 to 183	64-102	73	39-54	57	M80	320-0045	✓
Under 165	64-102	65	39-54	57	M80S	320-0052	✓
183+	64-102	83	39-54	57	M80L	N/A	
165 to 183	Under 64	73	35-39	44	M76	333-7581	✓
Under 165	Under 64	65	35-39	44	M76S	333-7599	✓
183+	Under 64	83	35-39	44	M76L	N/A	
165 to 183	102-140	73	52-57	59	M180	320-0060	✓
Under 165	102-140	65	52-57	59	M180S	320-0078	✓
183+	102-140	83	52-57	59	M180L	N/A	
152 to 183	140 +	73	65-90	94	M200HL	N/A	✓



## Adult Full Leg

For casts or dressings covering the whole leg.

- Seals at the top of the leg
- Suitable for use with diabetic ulcers
- Suitable for the bath or shower
- NHS SBS framework options



Height (cm)				rence (cm)	Model	NHS SBS
,	112.8.12 (1.8/		Limb			
165 to 183	64-102	96	52-65	70	M100	✓
Under 165	64-102	85	52-65	70	M100S	
183+	64-102	103	52-65	70	M100L	
Under 183	Under 64	96	39-54	57	M80LL	✓
152 to 183	102+	96	65-90	94	M200FI	✓



## Adult Knee

For light dressings on the knee.

- Seals above and below the knee
- For shower use only



				ence (cm) Dressing	Model
165 to 183	64-102	38	39-54	57	M86

## Adult Orthopaedic Boot

For orthopaedic boots. For the seal to be large enough to go over the boot, please select the appropriate size Adult Full Leg LimbO model.

- Seals at the top of the thigh
- Suitable for the bath or shower



## **Adult Foot**

For light dressings below the ankle.

- Seals above the ankle bone
- Suitable for use with diabetic ulcers
- For shower use only
- Available on FP10 prescription
- NHS SBS framework options



Shoe Size	Circumfer Limb	rence (cm) Dressing	Model	PIP Code	NHS SBS
Max shoe size 11	20-25	33	M20	416-6476	✓
Max shoe Size 11	23-34	40	M25	416-6484	✓

# Eliminate the need for replacing wet bandages or ruined plaster casts



## Child Full Leg

For casts or dressings covering the leg.

- Seals at the top of the thigh
- Suitable for the bath or shower



Age	Weight (kg)	Height	Length (cm)	Circumfer Limb		Model
3 yrs			53	25-34	40	FL3
4-5 yrs	Size is tv	pical for age,	58	25-34	40	FL45
6-7 yrs	select a	s you would	61	35-39	44	FL67
8-10 yrs	for	clothing	73	39-54	57	FL810
11-13 yrs			83	39-54	57	FL1113



## Child Half Leg

For casts or dressings below the knee, on older children.

- Seals just above the knee
- Suitable for the bath or shower



Age						Model
8-10 yrs	Size is typical	for age, select	58	25-34	40	BK810
11-13 yrs	as you would	d for clothing	61	35-39	44	BK1113



## Child Full Arm

For casts or dressings on the arm.

- Seals at the top of the arm
- Aqua blue junior style
- Suitable for the bath or shower



Age	Weight (kg)	Height	:	Length (cm)	Circumfe Limb	rence (cm) Cast	Model
2-3 yrs				31	16-22	29	FA23
4-5 yrs	Size is typ	ical for age,		43	16-22	29	FA45
6-7 yrs	select as	you would		49	22-25	33	FA67
8-10 yrs	for cl	othing		58	22-25	33	FA810
11-13 yrs				58	25-34	40	FA1113



## Child Half Arm

For casts or dressings below the elbow, on older children.

- Seals just above the elbow
- Slim and standard arm options
- Aqua blue junior style
- Suitable for the bath or shower



Age	Weight (kg)	Height	Length (cm)	Circumferen Limb	ice (cm) Cast	Model
8-10 yrs-Slim			36	16-22	29	BE810-40
8-10 yrs-Standard	Size is typical for		36	22-25	33	BE810-50
11-13 yrs-Slim	select as you v for clothin		43	16-22	29	BE1113-40
11-13 yrs-Standard		.0	43	22-25	33	BE1113-50

# Accessories and protection for casts, dressings and PICC Lines

## Making life more comfortable



## Arm Sleeve for Cast / Dressing / PICC Line

Brighten up your arm cast with a soft lycra cover.

- Helps avoid snagging
- · Adult and child sizes
- Elbow / PICC sleeve options smart and discreet









## Adult Leg Sleeve for Cast / Dressing

Keep your leg cast clean with a stylish lycra cover.

- Helps avoid snagging
- Available colours



## Toe Cozy

Keep toes clean and warm with a cosy fleece cover.

- Adjustable Velcro® straps
- Can be used with dressings too







## Oarsome Grips

Extra comfort for crutch handles.

- Helps avoid blisters
- Sold as a pair
- Available colours





## Keep clean and dry whilst outside

Adult Foot Cover (Sandal not included)

Keep your cast / dressing clean and dry whilst outside.

- Breathable waterproof fabric
- Adjustable Velcro® straps
- Strong and durable



Adult Mitt / Child Mitt

Keep your cast / dressing clean and dry whilst outside.

- Breathable waterproof fabric
- Adjustable Velcro<sup>®</sup> straps
- Strong and durable



Adult Wrist Cover

Keep your cast / dressing clean and dry whilst outside with exposed fingers for a better grip.

- Breathable waterproof fabric
- Adjustable Velcro® straps
- Left and right hand options





Rugged outdoor sock for child's lower leg.

- Perfect for the garden or beach
- Waterproof and breathable



# Diabetic Sock Collection - Socks that do not constrict blood flow on the lower leg

#### **Prosox**

- Stay up without applying pressure
- Welt-free
- Available in cotton or wool
- Anti-slip option available

## **Gentle Grip Socks**

- Non-binding
- Available in black, grey or six vibrant colours as shown
- Extra wide oedema option available





## Sock Aid

- Pull on socks with minimum stretching
- Ribbed thumb grips
- Lightweight



## Walk-Ease Walking Stick

- Adjustable length
- Handle designed for extra comfort
- Flexible ferrule to maximise contact with the ground
- Three styles

## Stick and Stay (Suction Grab Rail)

- Easy to fit
- Ideal for travel
- Suitable for ceramic or glass
- Safety indicator



## Helping you whilst you help patients

The bodyblock limb support provides a hygienic multidisciplinary, multipurpose support block for clinical environments.



With four sides creating different heights, the

bodyblock is a versatile tool for delivering stable elevation during treatment. Simply flip it for a different height and reduce bending to treat patients. The bodyblock

can be used to support the head, legs, heels or arms and can also be used as knee spacer. The bodyblock is lightweight, strong and easy to sterilise.



#### NO OLUBBI E GLIARANTEE

We want users to be totally happy with their LimbO and the service they receive from us. The LimbO is guaranteed to last twelve weeks when used correctly. If the LimbO proves unsuitable for the intended application or fails to last twelve weeks, a refund will be given or replacement will be provided. This does not affect statutory rights.

Please note: The LimbO comes with full instructions and safety information The LimbO is not intended for use on open or undressed wounds and should not be used in the sea or swimming pools.

Prices correct at September 2021

12 13


## Adult (for quick size selection use height and weight. Use length and circumference points shown on diagram for extra guidance)

	Limb Coverage	Full Product Description	Height (cm)	Weight (kg)	Limbo Length (cm)	Limb Circumference (cm)	Cast Circumference (cm)	Model	PIP Code	NHS SBS	Diabetic Ulcers	PICC Lines	Page	Conversion Tables	
Ī	Full	Adult Full Arm-Slim	152+	51-70	69	25-29	40 (maximum)	M60L		✓			4	Stone	Kilos
	Full	Adult Full Arm-Medium	152+	70-108	69	30-39	44 (maximum)	M70		✓			4	8	51
	Full	Adult Full Arm-Large	152+	108+	69	39-54	57 (maximum)	M87		✓			4	8.5	54
	Half	Adult Half Arm-Slimmer Adult	152+	Under 51	49	22-25	33 (maximum)	M50A					4	9.5	57 60
	Half	Adult Half Arm-Slim to Medium	152+	51-70	53	25-29	40 (maximum)	M60		✓			4	10	64
	Half	Adult Half Arm-Medium to Large	152+	70+	53	30-39	44 (maximum)	M67		✓			4	10.5	67
	Elbow	Adult Elbow-Extra Slim	152+	Under 51	26	22-25	33 (maximum)	M45				•	5	11	70
	Elbow	Adult Elbow-Slim to Medium	152+	51-70	26	25-29	40 (maximum)	M65	386-4303	✓		•	5	11.5 12	73 76
	Elbow	Adult Elbow-Medium to Large	152+	70-108	26	30-39	44 (maximum)	M75	386-4311	✓		•	5	12.5	79
	Elbow	Adult Elbow-Extra Large	152+	108+	26	39-54	57 (maximum)	M85				•	5	13	83
	Full	Adult Hand	Slim to M	edium Build	30	22-25	Bandage Only	MITT					5	13.5	86
	Full	Adult Full Leg	165 to 183	64-102	96	52-65	70 (maximum)	M100		✓			6	14 14.5	89 92
	Full	Adult Full Leg-Short	Under 165	64-102	85	52-65	70 (maximum)	M100S					6	15	95
	Full	Adult Full Leg-Long	183+	64-102	103	52-65	70 (maximum)	M100L					6	15.5	98
	Full	Adult Full Leg-Slim	Under 183	Under 64	96	39-54	57 (maximum)	M80LL		✓			6	16	102
	Full	Adult Full Leg-Extra Heavy Leg	152 to 183	102+	96	65-90	94 (maximum)	M200FL		✓			6	16.5 17	105 108
	Half	Adult Orthopaedic Boot	For any Orth	opaedic Boot please	e select a full leg LimbO	)							6	17.5	111
	Half	Adult Half Leg	165 to 183	64-102	73	39-54	57 (maximum)	M80	320-0045	✓	•		6	18	114
	Half	Adult Half Leg-Short	Under 165	64-102	65	39-54	57 (maximum)	M80S	320-0052	✓	•		6	18.5	118
П	Half	Adult Half Leg-Long	183+	64-102	83	39-54	57 (maximum)	M80L			•		6	19 19.5	121 124
1	Half	Adult Half Leg-Slim	165 to 183	Under 64	73	35-39	44 (maximum)	M76	333-7581	✓	•		6	20	127
	Half	Adult Half Leg-Slim-Short	Under 165	Under 64	65	35-39	44 (maximum)	M76S	333-7599	✓	•		6	Feet	
	Half	Adult Half Leg-Slim-Long	183+	Under 64	83	35-39	44 (maximum)	M76L			•		6	4	<b>cm</b> 122
	Half	Adult Half Leg-Large	165 to 183	102-140	73	52-57	59 (maximum)	M180	320-0060	✓	•		6	4,6	137
	Half	Adult Half Leg-Large-Short	Under 165	102-140	65	52-57	59 (maximum)	M180S	320-0078	✓	•		6	5	152
	Half	Adult Half Leg-Large-Long	183+	102-140	83	52-57	59 (maximum)	M180L			•		6	5,6	168
	Half	Adult Half Leg-Extra Heavy Leg	152 to 183	140 +	73	65-90	94 (maximum)	M200HL		✓	•		6	6,6	183 198
	Knee	Adult Knee	165 to 183	64-102	38	39-54	57 (maximum)	M86					7	7	214
)	Full	Adult Foot-Slim to Medium	Max shoe size 11			20-25	33 (maximum)	M20	416-6476	✓	•		7		
	Full	Adult Foot-Medium to Large		Max shoe Size 1	1	23-34	40 (maximum)	M25	416-6484	✓	•		7		

## Child (for quick size selection use clothing size. Use length and circumference points shown on diagram for extra guidance)

Cillic	√ (for quick size selection use clothing size. Use)	e lengui anu	a circumierer	ce points sno	iwn on diagram for extra	guidance)					
Full	Child Full Arm-2-3 yrs				31	16-22	29 (maximum)	FA23			9
Full	Child Full Arm-4-5 yrs				43	16-22	29 (maximum)	FA45			9
Full	Child Full Arm-6-7 yrs				49	22-25	33 (maximum)	FA67			9
Full	Child Full Arm-8-10 yrs		Size is typical for		58	22-25	33 (maximum)	FA810			9
Full	Child Full Arm-11-13 yrs		age, sele as you woul		58	25-34	40 (maximum)	FA1113			9
Half	Child Half Arm-8-10 yrs-Slim		clothing		36	16-22	29 (maximum)	BE810-40			9
Half	Child Half Arm-8-10 yrs-Standard				36	22-25	33 (maximum)	BE810-50			9
Half	Child Half Arm-11-13 yrs-Slim				43	16-22	29 (maximum)	BE1113-40			9
Half	Child Half Arm-11-13 yrs-Standard				43	22-25	33 (maximum)	BE1113-50			9
<b>/</b> Elbow	Child Elbow	Under 152		Under 51	12	22-25	33 (maximum)	M45S		•	5
Full	Child Full Leg-3 yrs				53	25-34	40 (maximum)	FL3			8
Full	Child Full Leg-4-5 yrs				58	25-34	40 (maximum)	FL45			8
Full	Child Full Leg-6-7 yrs		Size is typic		61	35-39	44 (maximum)	FL67			8
Full	Child Full Leg-8-10 yrs		age, sele		73	39-54	57 (maximum)	FL810			8
Full	Child Full Leg-11-13 yrs		as you would for clothing		83	39-54	57 (maximum)	FL1113			8
Half	Child Half Leg-8-10 yrs				58	25-34	40 (maximum)	BK810			8
Half	Child Half Leg-11-13 yrs				61	35-39	44 (maximum)	BK1113			8

## Measuring Guide



# Here to help you, your team and your patients









